



## Genomtec® SARS-CoV-2 RT-LAMP/N CE-IVD Laboratory Kit

# Instructions for use

Real-Time Reverse Transcription Loop-Mediated Isothermal Amplification Test for qualitative detection of nucleic acid from SARS-CoV-2 virus.



Genomtec SA, ul. Stabłowicka 147, 54-066, Wrocław, Poland

Tel: +48 793 440 931

---

The customer is responsible for compliance with regulatory requirements that pertain to their procedures and uses of the instrument.

The information in this guide is subject to change without notice.

Please ensure that you are using a current version of the IFU. Refer to <http://genomtec.com/support> for latest version.

Document revision	Date	Description
1	27 May 2020	Initial release

---

## Table of contents

1. Definitions	5
2. References	6
3. Genomtec® SARS CoV-2 RT-LAMP/N laboratory kit general product information	7
3.1. Intended purpose	7
3.2. Summary and product description	7
3.3. Principles of the procedure	7
3.4. Reagents and materials	9
3.4.1. Materials and reagents provided	9
3.4.2. Materials required but not provided	9
4. Precautions	11
5. Storage, shipping and stability	12
6. Quality control	13
6.1. Assay controls	13
6.2. Specimen collection, handling, transport and storage	13
7. Operating instructions	14
7.1. Before starting the test	14
7.1.1. Preparation step 1	14
7.1.2. Preparation step 2	14
7.1.3. Preparation step 3	14
7.1.4. Preparation step 4	14
7.2. Performing the test	14
7.2.1. Test step 1	14
7.2.2. Test step 2	15
7.2.3. Test step 3	15
7.2.4. Test step 4	15
7.2.5. Test step 5	15
8. Result interpretation	17
8.1. Interpreting patients' results	17
8.2. Limitations	17

---

9. Performance characteristics	19
9.1. Analytical sensitivity (Limit of Detection)	19
9.2. Analytical reactivity (Inclusivity)	19
9.3. Analytical specificity (Cross Reactivity)	19
9.4. Clinical evaluation	19
9.4.1. Conclusion:	21
10. Symbols	22
11. Ordering and contact information	23

---

## 1. Definitions

Abbreviation	Definition
SARS-CoV-2	Severe Acute Respiratory Syndrome
LAMP	Loop-Mediated Isothermal Amplification
DNA	Deoxyribonucleic Acid
RNA	Ribonucleic Acid
RNase	Ribonuclease
DNase	Deoxyribonuclease
IFU	Instructions for Use
PI	Product Information
MERS	Middle East Respiratory Syndrome
NATT	Nucleic Acid Amplification Test
cDNA	Complimentary DNA
FAM	Fluorescein amidite
PCR	Polymerase Chain Reaction
LOD	Limit of detection
RT	Reverse transcription
BLAST	Basic Local Alignment Search Tool

## 2. References

Tsugunori Notomi, Hiroto Okayama, Harumi Masubuchi, Toshihiro Yonekawa, Keiko Watanabe, Nobuyuki Amino, and Tetsu Hase. Loop-mediated isothermal amplification of DNA. *Nucleic Acids Res.* 2000 Jun 15; 28(12): e63.

Hong TC, Mai QL, Cuong DV, et al. Development and evaluation of a novel loop-mediated isothermal amplification method for rapid detection of severe acute respiratory syndrome coronavirus. *J Clin Microbiol.* 2004;42(5):1956-1961. doi:10.1128/jcm.42.5.1956-1961.2004

### 3. Genomtec® SARS CoV-2 RT-LAMP/N laboratory kit general product information

#### 3.1. Intended purpose

The Genomtec® SARS CoV-2 RT-LAMP/N test is a CE-IVD Laboratory Kit containing controls and reagents intended for reverse transcription and amplification of nucleic acid in an isothermal reaction, specifically Loop-Mediated Isothermal Amplification (LAMP). It is a qualitative assay detecting specifically SARS-CoV-2 in throat swab and nasopharyngeal swab specimens from individuals suspected of COVID-19.

The results obtained will identify the presence of SARS-CoV-2 RNA in the sample. The positive result obtained with the diagnostic test should only be taken into consideration together with patient's clinical history and other diagnostic results while concluding the final infection status.

Similarly, negative results do not exclude entirely COVID-19 disease and should be accompanied by other diagnostic solution ruling out absence of SARS-CoV-2 in patients. The wholistic patient management should carefully consider spectrum of clinical symptoms, patient history, linked to available epidemiological data.

Testing with the Genomtec® SARS-CoV-2 RT-LAMP/N Laboratory Kit is intended for molecular diagnostic clinical laboratory use by qualified and trained clinical laboratory personnel.

- Note:- The test is not sterile, and does not require a sterile operating environment.
- Note: - The assay is not for self testing.

#### 3.2. Summary and product description

Genomtec® RT-LAMP/N Laboratory Kit contains assay enough to perform 50 reactions (including controls and assay mixes) required for the RT-LAMP detection of RNA from the SARS-CoV-2. Specifically, Genomtec® RT-LAMP/N Laboratory Kit confirms presence of N gene encoding Nucleocapsid spike protein in the analysed human sample.

Genomtec® SARS-CoV-2 RT-LAMP/N CE-IVD Laboratory Kit includes the following reagents:

<b>Mastermix</b>	Vial containing Genomtec® SARS-CoV-2 Detection Mix (50 reactions) with primers recognising specific fragment of N gene SARS-CoV-2.
<b>Control   I</b>	Vial containing Genomtec® SARS-CoV-2 Inhibition Control Mix (50 reactions) with primers recognising specific fragment of the human genome controlling appropriateness of biological sample collection and RNA purification.
<b>Control   +</b>	Vial containing Genomtec® SARS-CoV-2 Positive Control in the form of a synthetic SARS CoV-2 cDNA.
<b>D. Water</b>	Vial containing DNase/RNase-Free Distilled Water.

The kit also includes the product information insert (PI00A) which provides the instructions and the download link for the Instructions for Use (this document).


#### 3.3. Principles of the procedure

Genomtec® RT-LAMP/N Laboratory Kit targets distinct genomic regions of SARS-CoV-2 N gene that is conserved and are present as the N gene is present in coronaviruses (e.g., SARS, MERS, bat SARS-like coronavirus, etc.).

Genomtec's recommendation for use Product line Genomtec® SARS-CoV-2 RT-LAMP/N CE-IVD kit are as follows:

- Where new or suspected cases of COVID-19 disease arise in the population utilise Genomtec® SARS-CoV-2 RT-LAMP/N CE-IVD Kit to aid in diagnosis of SARS-CoV-2. Positive result for SARS CoV-2 (detected gene N) and all assay controls indicate active infection.
- False negative result (all controls negative; please see Section 8, Results Interpretation) requires patient's resampling (recommended) subject to nucleic acid purification and test repeat using either Genomtec® SARS-CoV-2 RT-LAMP/N CE-IVD kit or other manufacturer's genetic assay.
- If the second diagnostic round confirms presence of both SARS-CoV-2 RNA and the manufacturers controls, the result is positive. If results are discordant then manufacturer's confirmatory protocol must be executed. It may include if appropriate, sequencing of the virus from the original specimen or an amplicon generated from an appropriate NAAT assay.

Please refer to the schematic diagram 'COVID-19 diagnostic pathway'





Genomtec recognises difficulties in recognition of a proper epidemiological situation and way of SARS-CoV-2 virus transmission in population, therefore the User is advised to follow the latest guidance provided by World Health Organisation (WHO).

The appropriate decision on epidemic status and virus transmission should be provided by each Country's appropriate Healthcare and / or Epidemiological Agencies. Genomtec will not provide any advice on epidemiological status of any country or geographical area.

The diagnostic workflow consists of:

- Sample collection and nucleic acid (RNA) purification.
- Reverse transcription of the purified nucleic acid and simultaneous cDNA amplification using the Genomtec® SARS-CoV-2 RT-LAMP/N CE-IVD kit.
- The RT-LAMP process occurs at constant temperature (isothermal condition), where set of primers (five primers) recognises seven specific genomic sequences encompassing targeted SARS-CoV-2 gene N. In the first stage Reverse Transcription is being performed generating cDNA. Simultaneously with RT in the same constant temperature all five primers anneal and amplify the cDNA template. In presence of targeted genomic fragment increasing cDNA concentration prompts fluorescent dye binding, simultaneously increasing the fluorescence intensity. As LAMP reaction is not exponential, increasing-cycle based fluorescence might not be expected (fluorescence may occur suddenly, without the correlation between higher cDNA initial concentration and stronger and / or quicker fluorescent signal achieved).
- Fluorescence is monitored by the RealTime-PCR instruments (equipped with FAM channel detector) and the readout is provided that next undergoes analysis (see Section 8, Interpretation of Results).

### 3.4. Reagents and materials

#### 3.4.1. Materials and reagents provided

<b>Mastermix</b>	Vial containing Genomtec® SARS-CoV-2 Detection Mix (50 reactions) with primers recognising specific fragment of N gene SARS-CoV-2.
<b>Control  </b>	Vial containing Genomtec® SARS-CoV-2 Inhibition Control Mix (50 reactions) with primers recognising specific fragment of the human genome controlling appropriateness of biological sample collection and RNA purification.
<b>Control +</b>	Vial containing Genomtec® SARS-CoV-2 Positive Control in the form of a synthetic SARS CoV-2 cDNA.
<b>D. Water</b>	Vial containing DNase/RNase-Free Distilled Water.

#### 3.4.2. Materials required but not provided

- Thermal Cycler with FAM channel for fluorescence detection and 96 well plate holder or block suitable for strips of PCR-tubes with optically clear caps, maintained and calibrated according the the manufacturers instructions.
- Laboratory freezers: - 30°C to -10°C
- Laboratory microcentrifuge (for 48 microtubes size 1.5 to 2 mLs) and a centrifuge with a rotor for microplates (speed: up to 17,500 rpm)

- Laboratory mixer, vortex or equivalent
- Single and / or multichannel adjustable pipettes working in volume range:
  - 0,5-10µl
  - 10-100µl
  - 100-1000µl
- Cooling block or ice
- PCR 96-Well Reaction Plate or strips of PCR-tubes with optically clear caps
- Optical Adhesive Film
- Sterilise aerosol barrier (filtered) pipette tips

---

## 4. Precautions

Testing with the Genomtec® SARS-CoV-2 RT-LAMP/N CE-IVD Laboratory Kit is intended for use by qualified and trained clinical laboratory personnel specifically instructed and trained in the techniques of RealTime-PCR and / or LAMP and in vitro diagnostic procedures. Use separate areas for the preparation of patient samples and controls to prevent cross-contamination.

Samples and reagents must be handled under a PCR laminar airflow hood or biological safety cabinet. PCR workstation should be sterilised with UV light for minimum of 15 minutes before use.

- All specimens should always be treated as potentially infectious and/or biohazardous in accordance with safe laboratory procedures.
- Use personal protective equipment (PPE) according to local guidelines for the handling of potentially infectious samples.
- Always use sterile, nuclease-free pipette tips with aerosol barriers (filtered).
- Do not eat, drink or smoke in the working area.
- Manufacturer does not provide warranty if any modifications to assay reagents, assay protocol, or instrumentation were made, and these are in violation of the In Vitro Diagnostic Directive 98/79/EC.
- Do not use the kit after the expiry date.
- Never open the plate or PCR-tubes after the analysis process as it might cause contamination by the amplicons present in positive reactions.
- Dispose of waste in compliance with the local biohazard regulations. Check safety procedures set by your institution for working with chemicals and handling biological specimens.
- Avoid sample or reagent contact with the skin, eyes, and mucous membranes. If skin, eyes, or mucous membranes come into contact, rinse immediately with water and seek medical advice immediately.
- In case of sample or reagents coming in contact with skin, eyes or mucous membranes, or if swallowed, immediately follow the laboratory post-exposure protocol.
- Clean and disinfect all reagents and / or sample spillage with disinfectant such as 0.5% sodium hypochlorite, or other suitable disinfectant
- Safety Data Sheets are available upon request with Genomtec or Authorised Distributor.
- Laboratories may be required to report all positive results to the appropriate local Health Authorities.
- Positive results in this test indicate presence of SARS-CoV-2 RNA in a patient sample.
- Reagents must be stored and handled as specified in Table 1.
- The quality of RNA preparation (purification) may influence the quality of the RT-LAMP reaction, therefore it is recommended for laboratories to use previously validated purification method.

## 5. Storage, shipping and stability

IMPORTANT! POSITIVE CONTROL SHOULD BE ALIQUOTED INTO SMALLER VOLUMES TO PREVENT ITS DEGRADATION AND PROTECT AGAINST MULTIPLE THAW-FREEZING CYCLES.

Detection mix and inhibition control mix are stable for at least three freeze-thaw cycles.

Table 1.

Reagent	Quantity	Volume	Storage limits	Shipping limits	Shelf life
Genomtec® SARS-CoV-2 Detection Mix	1 vial	750 µl	-22°C to -15°C	Dry ice	Three months
Genomtec® SARS-CoV-2 Inhibition Control Mix	1 vial	750 µl	-22°C to -15°C	Dry ice	Three months
Genomtec® SARS-CoV-2 Positive Control	1 vial	40 µl	-22°C to -15°C	Dry ice	Three months
DNase/RNase-Free Distilled Water	1 vial	1000 µl	-22°C to -15°C	Dry ice	Three months

## 6. Quality control

### 6.1. Assay controls

Positive, negative test controls as well as the inhibition control must be included to accurately interpret patient test results. Including the inhibition control minimises occurrence of potential false negatives.

Include the below Controls:

Type of Control	Contents and targets	Function
Positive Control	Synthetic SARS-CoV-2 cDNA with assay primers directed against specific genomic sequence of targeted N gene	Ensure the proper reaction conditions as well as stability of the Reaction mix
Negative Control	Reaction mix with primers recognising specific sequences on targeted gene N of SARS-CoV-2 combined with DNase/ RNase-Free Distilled Water	Ensures lack of cross-contamination arising from assay set-up
Inhibition Control	Inhibition control mix with added patient's extracted RNA sample	Controls possible inhibition of the amplification and appropriate sample collection procedure (e.g. by throat swab taken from human) and RNA purification efficiency

### 6.2. Specimen collection, handling, transport and storage

Patient samples must be collected according to appropriate laboratory guidelines. These include throat swab and nasopharyngeal swab specimens.

Treat all samples and controls as if they are capable of transmitting infectious pathogen.

Samples should be transported on Universal Transport Medium, or other recommended equivalent, and stored at -20°C / -70°C. Samples RNA must be shipped on dry ice.

- Note: nasopharyngeal aspirate, and bronchoalveolar lavage (BAL) specimens have not been validated for use

---

## 7. Operating instructions

### 7.1. Before starting the test

#### IMPORTANT

- If working with large number of samples, to minimise degradation of RNA analytes keep the plate / PCR microtubes on ice / in cooling block until it is loaded into the RealTime-PCR instrument.
- Use thermocycler to run the plate immediately after preparation. Failure to do so could result in degraded RNA samples.
- Prevent contamination implementing separate areas for RNA purification and reaction amplification; prepare reagents in a PCR workstation (with dual decontamination action by UV) and use separate pipette for controls and samples, and always use aerosol barrier pipette tips
- Maintain an RNase-free environment.
- Protect kit components (particularly Assay Mix and Internal Control Mix) from light.
- Each sample requires one Inhibition Control to be run and positive and negative controls are required to be run for each assay.( see “Quality Control” in Section 6)
- After thawing please place all reagents and assay components onto the cooling block or ice to preserve their potency.

#### 7.1.1. Preparation step 1

The first step in the diagnostic workflow is sample collection (from throat swab or nasopharyngeal swab specimens) and sample RNA purification.

#### 7.1.2. Preparation step 2

The quality of RNA purification may influence the quality of the RT-LAMP reaction; therefore, it is recommended for laboratories to use previously validated and commercially available purification method.

#### 7.1.3. Preparation step 3

All provided and necessary reagents and assay components should be defrosted at  $\leq 4^{\circ}\text{C}$ , followed by gentle vortexing (mixing) and brief centrifugation (to collect evaporate and liquid from the cap).

#### 7.1.4. Preparation step 4

Set-up the reaction according to Point 7.2 Performing the test

### 7.2. Performing the test

#### 7.2.1. Test step 1

Prepare all assay Controls and samples according to the Table 2. Each sample requires simultaneous run with assay Controls (Positive, Negative, Inhibition). Each run of the process requires additional reaction for positive control and negative control to be added. Pipette to designated wells Detection Mix, Inhibition Control Mix and DNase/RNase-Free Distilled Water, and then add either the Sample or the required Control. The final volume in each well is 20  $\mu\text{l}$ .

Table 2. Reaction plate set-up.

Reagent	Analyte	Inhibition Control	Positive Control	Negative Control
Genomtec® SARS-CoV-2 Detection Mix	15 µl	-	15 µl	15 µl
Genomtec® SARS-CoV-2 Inhibition Control Mix	-	15µl	-	-
Sample RNA	5 µl	5 µl	-	-
Genomtec® SARS-CoV-2 Positive Control	-	-	5 µl	-
DNase/RNase-Free Distilled Water	-	-	-	5 µl
Total Volume	20 µl	20 µl	20 µl	20 µl

**7.2.2. Test step 2**

After each addition, pipette up and down ensuring proper mixing.

**7.2.3. Test step 3**

Seal the plate with Optical Adhesive Film, then centrifuge briefly to collect the liquid at the bottom of the reaction plate.


**7.2.4. Test step 4**

Place in the RealTime-PCR instrument that is configured as follows:

Step	Temperature [°C]	Time [Sec.]	Cycles/ repeats
Amplification 1	62	30	30
Amplification 2	62	30	

**7.2.5. Test step 5**

The RealTime-PCR instrument must be able to operate on 20 µl total volume in an individual well of a multi-well PCR plate or individual tube in a strip of PCR-tubes with optically clear caps and must be equipped with optics and filters allowing fluorescence reading at FAM/SYBR/Green channel. Use provided Instrument's Software to set up the run protocol. Below you can find an example of the S-curves achieved for a positive and negative samples, and also the threshold.





## 8. Result interpretation

### 8.1. Interpreting patients' results

To interpret results of the assay please follow the guidance presented in Table 3.

Table 3. Result interpretation for patient samples

The Analyte	Inhibition Control	Positive Control	Negative Control	Expected Result
+	+	+	-	POSITIVE
-	+	+	-	NEGATIVE
+	+	+	+	FALSE POSITIVE
-	-	-	-	FALSE NEGATIVE

#### ACTIONS:

- **POSITIVE** - report results to the appropriate Local Health Agency / Provider, as applicable.
- **NEGATIVE** - report results to the healthcare provider; if patient in clinical review symptomatic, consider other diagnostic tests for other pathogens.
- **FALSE POSITIVE** - Repeat test on used RNA sample and if the repeat result comes inconclusive, consider collecting a fresh specimen for retesting.
- **FALSE NEGATIVE** - Repeat test (recommended fresh resampling) using a genetic diagnostic kit. If the repeat result remains discordant, additional confirmatory genetic testing should be conducted, particularly in presence of clinical symptoms.

### 8.2. Limitations

Genomtec® SARS-CoV-2 RT-LAMP/N Laboratory Kit is intended for molecular diagnostic clinical laboratory use by qualified and trained clinical laboratory personnel. The Laboratory should have implemented Quality System and work in accordance to GLP, and in compliance with the guidelines presented in this Document in order to prevent cross-contamination of RNA clinical samples and other components of the Kit.

- Genomtec® SARS CoV-2 RT-LAMP/N Laboratory Kit was validated on throat and nasopharyngeal swabs for SARS-CoV-2 detection.
- All specimen collection, shipment and storage must be performed according to Section 5 of this Document and country specific guidelines for biological material handling and storage. All reagents and assay components must be stored in accordance to conditions described in Table 1. Failure to comply with the guideline may negatively affect the diagnostic procedure, providing false results.
- It is mandatory to implement previously laboratory-validated RNA purification method to provide the highest quality of sample RNA.
- The FALSE NEGATIVE results may be indicative of:
  - Unsuccessful biological material collection
  - RNA material degradation during storage / transport (lack of guideline compliance)
  - Inefficient RNA purification

- 
- Presence of RT-LAMP inhibitors in the reaction
  - Mutation in targeted SARS-CoV-2 amplicon (although if point mutation, highly unlikely to cause primer mismatch due to their mismatch-recognition ability)
  - Lack of compliance and execution of the diagnostic stages according to the enclosed Document
  - The FALSE POSITIVE results may be indicative of:
    - Cross-contamination of the assay components and / or samples (RNA) during reaction set-up
    - Mix-up of samples
    - Extrinsic RNA contamination during set-up procedure
  - The negative result obtained with the diagnostic test (Genomtec® SARS-CoV-2 RT-LAMP/N CE-IVD Kit) does not disqualify presence of SARS-CoV-2 and any diagnostic recommendation should not be made solely on its basis, instead it is recommended to consider other diagnostic results and patient's clinical history while concluding the final infection status.
  - Laboratories may have to report all positive results to the appropriate Competent Health Authorities.

## 9. Performance characteristics

### 9.1. Analytical sensitivity (Limit of Detection)

The LOD protocol defined the lowest number of copies of SARS-CoV-2 gene N that can be detected by the Genomtec® SARS-CoV-2 RT-LAMP/N CE-IVD Laboratory Kit.

Limit of detection of the assay was determined by performing reactions for series of dilutions of SARS-CoV-2 synthetic full genome RNA control and calculated in probit analysis. The LOD of the SARS-CoV-2 N gene assay is: 147 (95% CI: 112 - 271) gene copies / reaction.

### 9.2. Analytical reactivity (Inclusivity)

The assay primers were mapped to SARS-CoV-2 reference genome sequence NC\_045512.2:28274-29533 Severe acute respiratory syndrome coronavirus 2 isolate Wuhan-Hu-1, complete genome. LAMP primers (five) sequences (seven) for SARS-CoV-2 N gene assay showed 100% homology to all SARS-CoV-2 isolate analysed. Primers are designed to amplify conservative sequence of N gene of SARS-CoV-2.

In order to confirm coverage by the primers an alignment of 45 sequences of the N gene from the whole genome sequencing analysis of SARS-CoV-2 has been performed. 100% inclusivity was confirmed for 6 out of 7 regions of the amplified fragment of N gene (mutation site was detected in the 5' end of one primer, which due to the nature of the amplification technology does not influence the reaction efficiency). Five primers were used further in the final assay.

### 9.3. Analytical specificity (Cross Reactivity)

Specificity of Genomtec® SARS-CoV-2 RT-LAMP/N CE-IVD Laboratory Kit was confirmed in a study where assay mix (containing assay primers) was spiked with the below potentially cross-reacting pathogens and the sample was subjected to RT-LAMP run. Furthermore, additional in-silico analysis of designed primers utilising BLAST alignment tool was conducted and show no possible full sequence similarity

The twenty four (24) pathogens utilised in this study included:

Mycoplasma genitalium	Escherichia coli
Streptococcus pyogenes	Candida albicans
Enterococcus faecalis	Mycoplasma pneumoniae
Moraxella catarrhalis	Klebsiella pneumoniae
Legionella pneumophila	Staphylococcus aureus methicilin sensitive (MSSA)
Enterococcus faecium	Acinetobacter baumannii
Mycoplasma hominis	Ureoplasma urealyticum
Haemophilus influenzae	Human genomic DNA
Bordetella pertussis	Staphylococcus aureus methicilin resistant (MRSA)
Pseudomonas aeruginosa	Listeria monocytogenes
Haemophilus ducreyi	Campylobacter jejuni
Chlamydia pneumoniae	Mobiluncus mulieris

None of the above listed pathogens had any effect on RT-LAMP assay performance and cross-reacted.

### 9.4. Clinical evaluation

A clinical study was performed to evaluate the performance of the Genomtec® SARS-CoV-2 RT-LAMP/N CE-IVD Laboratory Kit on mixture of throat and nasopharyngeal swabs (NP) obtained from patients suspected of contracting / developing COVID-19 disease on territory of Republic of Poland. The study samples used in the investigation included leftover specimens collected for routine clinical care or analysis that would otherwise had been discarded. The pre-released version of test was evaluated by two independent medical laboratories in Poland, and a total of 43 clinical specimens were tested in April and May 2020 during the course of the study.

The method of sample collection and transport medium used (type and manufacturer) complied with the general clinical laboratory regulations enforced by relevant Polish Healthcare Authority (Panstwowy Zakład Higieny), and also with WHO recommendations.

Patients that underwent swab collection were referred to such diagnostic procedure based on either clinical symptoms, suspicion of contracting the disease due to contact with infected individual(s) or travelled from abroad and were subjected quarantine. All patients were referred to conduct the genetic diagnostic test against SARS-CoV-2 by a healthcare professional.

There was one swab collected from each patient that had priority to undergo testing utilising the standard care molecular diagnostic test based on a real time qPCR technology. If the remaining leftover RNA sample contained enough genetic material to set-up an assay using Genomtec® SARS-CoV-2 RT-LAMP/N CE-IVD Laboratory Kit, it was processed and executed on a standard RealTime-PCR instrument.

The qRT-PCR CE-IVD diagnostic kits were used as reference method for RT-LAMP assay and the protocol was compliant with the manufacturer's instructions. Investigated samples included in the study were also processed with the pre-released version of Genomtec® SARS-CoV-2 RT-LAMP/N CE-IVD Laboratory Kit according to instruction enclosed in this document. The nucleic acid purification was carried with one of the commercially available kits. A total of 35 specimens were included in the analysis of clinical performance. In 8 cases the specimens that produced unresolved results by RT-qPCR were excluded from the analysis (lack of confirmation for the presence of one gene out of two included in the reference assay). Table 4 depicts the overall assay clinical performance vs standard practice RT-qPCR method whereas table 5 presents results by site.

Table 4.

		Reference standard practice RT-qPCR assay result		
		Positive	Negative	Unresolved
<b>Genomtec® SARS-CoV-2 RT-LAMP/N CE-IVD Laboratory Kit</b>	<b>Positive</b>	32	-	32
	<b>Negative</b>	-	3	3
	<b>Unresolved</b>	32	3	35
Sensitivity (SE)		32/32 = 100% (95% CI :100%-100%)		
Specificity (SP)		3/3 = 100% (95% CI :100%-100%)		
Positive Predictive Value (PPV)		32/32 = 100% (95% CI :100%-100%)		
Negative Predictive Value (NPP)		3/3 = 100% (95% CI :100%-100%)		















Table 5. Stratified by site.

Site	Samples (%)	Reference RT-qPCR		Genomtec® SARS-CoV-2 RT-LAMP/N CE-IVD Laboratory Kit		SE	SP	PPV	NPV
		Positive	Negative	Positive	Negative				
A	32 (91%)	29 (91%)	3 (9%)	29 (19%)	3 (9%)	100%	100%	100%	100%
B	3 (9%)	3 (100%)	0	3 (100%)	0	N/A	100%	N/A	100%

**9.4.1. Conclusion:**

It has been confirmed that Genomtec(R) SRAS-CoV-2 RT-LAMP/N test exhibits 100% sensitivity and specificity compared with a standard (two gene) laboratory RT-qPCR diagnostic test when detecting presence or absence of the SARS-CoV-2 virus in clinical specimens. Both, the PPV and NPP were also obtained at level 100%.

## 10. Symbols

Symbol (IEC 15223-1 2016)	Description
	Indicates the need for the user to consult the instructions for use.
	Indicates the manufacturer's batch code so that the batch or lot can be identified.
	Indicates the total number of IVD tests that can be performed with the IVD.
	Indicates the manufacturer's catalogue number so that the medical device can be identified.
	Indicates a medical device that is intended to be used as an <i>in vitro</i> diagnostic medical device.
	Indicates the date after which the medical device is not to be used.
	Indicates the medical device manufacturer, as defined in EU Directives 90/385/EEC, 93/42/EEC and 98/79/EC.
	Indicates the temperature limits to which the medical device can be safely exposed.
	Indicates a medical device that should not be used if the package has been damaged or opened.
	Indicates a medical device that needs protection from light sources.
	Indicates the Genomtec® SARS-CoV-2 Detection Mix.
	Indicates the Genomtec® SARS-CoV-2 Positive Control.
	Indicates the Genomtec® SARS-CoV-2 Inhibition Control Mix.
	Indicates DNase/RNase-Free Distilled Water.

---

## 11. Ordering and contact information

Product	Order number
Genomtec® SARS CoV-2 RT-LAMP/N Laboratory Kit - UK language	GA00AUK

For ordering

Polorto: <http://vitea.cool/products/>

For technical support

Genomtec: <http://www.genomtec.com>